DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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If an area in the field of view is out of focus:

- ◆ Check the objective is properly located on the luminous beam's path. If not, turn it until it is in an appropriate position.
- ◆ Check the sample is properly located on the stage, making sure it is held by the stage clips. Then, if it does not focus revise the coarse adjustment knob.

If the frontal lens of the objective touches the prepared slide when focusing:

- ◆ Check the prepared slide is not the other way round (being the slide over the cover glass) and placing it properly in case it is appropriate to do so.

If the frontal lens of the more powerful objective touches the prepared slide when passing from a low-magnification objective to a superior one:

- ◆ Check the prepared slide is not the other way round and check also the slide's thickness, which has to have a standard thickness of 0.17 mm.

9. PRACTICAL RECOMMENDATIONS

- ◆ As the head is rotary you should rather observe by the stage's frontal part, as it facilitates access to the mechanical control knobs.
- ◆ Handling the stage: First, move the sample in consecutive steps, from the front to the back; each fraction is observed to detail by means of a lateral movement. Thus, the sample can be moved in a uniform way and so explore all its area.

INSTRUCTIONS ON ENVIRONMENTAL PROTECTION

Once the microscope has finished its life cycle, please do not throw it into the common trash, take it to a recycling point for electrical devices. The microscope has no dangerous or toxic elements for humans but an inappropriate disposal damages the environment.

Materials can be recycled. Recycling materials you contribute to protect the environment. Please, contact with your local entity to consult collection points.

Never adjust it with the condenser's diaphragm. Both the condenser and the lamp should be well centred. Check that between the field diaphragm and the aperture diaphragm there is not any diffuser filter.

- ◆ Eyepieces are well fit and objectives well screwed.
- ◆ Eyepieces and objectives are perfectly cleaned and free of particles. To check cleaning of the whole optical system proceed as follows:
 - Make the eyepieces turn checking the little specks are moving; if so, clean them. Make the head turn. You should never disassemble the head, but you can clean it with delicacy by blowing the accessible surfaces of prisms with a plastic bulb.
 - Turn the objective and, in case the parasite images also turn, clean it with the help of a dry brush to remove the dust. Watch the front part of the microscope with a magnifying glass or an inverted eyepiece.
- ◆ Thickness of slides, cover glasses and means of assembly is not too much as to avoid focusing with medium or high magnification (there are standard sizes for both the slides and cover glasses). Slides and cover glasses should be clean. Check that the slide is properly located and there are not two cover glasses superimposed.
- ◆ The combination objective-eyepiece is appropriate; the eyepiece can be too powerful for the selected objective.
- ◆ The prepared slide is well done; check it with a "test" prepared slide.

8. TROUBLESHOOTING

Before sending the microscope back to Zuzi Technical Assistance Department, check the following:

If the lamp does not switch on:

- ◆ Check if it is properly connected to the right current intake.
- ◆ Check if the lamp is in good conditions.

If the field of vision is cut:

- ◆ Check the nosepiece is properly fit. To do that, turn it slightly back and forth until it fits. Check it is well centred and if not, send it to Zuzi Technical Assistance Department through your distributor.

If there is dust or dirt in the field of view:

- ◆ Check there is dust on the precondenser lens, upper lens of the condenser or the eyepiece, as well as the cleanliness of the sample. Once you find out where the dirt is, clean it as explained before.

Thank you for choosing this equipment. We sincerely wish that you enjoy your Zuzi education microscope serie 116. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Zuzi educational microscope serie 116, so it should be available for all the users of this equipment.
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ You should carefully handle the microscope avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the microscope to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the microscope.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt, please, contact your distributor (installation, starting up, operation). Also you can send your doubts or suggestions to Zuzi Technical Support (asistencia@auxilab.es).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the microscope warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the microscope warranty.
- ◆ Neither fuses (1 A) nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Zuzi Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

 **ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

INDEX OF LANGUAGES

Spanish..... 2-14
 English..... 15-27
 French..... 28-39

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT 16
 2. DESCRIPTION..... 17
 3. TECHNICAL SPECIFICATIONS 17
 4. INSTALLATION / SETTING UP..... 19
 5. MAINTENANCE AND CLEANING..... 24
 6. CHOOSING OBJECTIVES AND EYEPIECES..... 25
 7. CAUSES OF A DEFECTIVE IMAGE 25
 8. TROUBLESHOOTING 26
 9. RECOMMENDATIONS 27
 ANNEX I: CE CERTIFICATE 40

1. USES OF THE INSTRUMENT

Thanks to their modern and developed design and to the modular configuration, Zuzi microscopes serie 116 becomes essential and ease-to-use equipment for laboratories in any education field.

the optical parts when it is not used for a long period of time.
 ♦ Please keep the original packaging either to transport the equipment, when it is not being used for a long time or when you send it for an overhaul.

Cleaning

- ♦ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the microscope and produce an early ageing of the equipment.
- ♦ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.
- ♦ Lenses must not be disassembled by the user. Were there any dust or dirt to be cleaned, clean it with a natural horse hair brush or a smooth piece of cloth, fluff-free, dampened with a bit of xilol or toluene.
- ♦ To remove dust from lenses, blow them with a plastic bulb or clean them with a natural horsehair brush.
- ♦ When cleaning mechanical or plastic parts do not use abrasive solvents that could damage either paint or finish. Always use neutral detergents.

 **ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

6. CHOOSING OBJECTIVES AND EYEPIECES

Zuzi microscopes serie 116 are supplied with WF10x eyepieces and with 4x, 10x, 40x and 100x achromatic objectives, but can be used with other eyepieces and objectives. To know our offer of eyepieces and objectives look up our catalogue or visit our web page (www.auxilab.com).

The observed image loses surface and sharpness as magnification increases. The mentioned increase should be done by changing objectives and putting each time a more powerful one and not by changing eyepieces to a higher magnification, as eyepieces only magnify the image obtained with the objective. Thus, the more magnification the eyepiece has, the higher will be the loss of sharpness, clarity and surface of the resulting image.

For routine observations use the eyepiece with lesser magnification with more powerful objectives. The eyepiece with the higher magnification should be kept back for particular occasions, bearing in mind that it decreases definition and does not increase resolution.

7. CAUSES OF A DEFECTIVE IMAGE

- In case of a defective image check the following:
- ♦ Illumination is well done and luminous intensity is neither excessive nor too weak.

Safety

- ◆ The microscope must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ Place the microscope on an horizontal, flat and stable table, making a free space for working.
- ◆ Do not place the microscope near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ When you are not using the microscope for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the microscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this microscope it is essential to follow the processes of use.


Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The microscope incorporate a supply adaptor, this adaptor should be at sight because in case of emergency you have to switch off.
- ◆ Never try to repair the microscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Zuzi Technical Assistance Department.
- ◆ Please use the plastic cover provided to keep the microscope away from dust lying on

2. DESCRIPTION

- | | |
|---------------------------|-------------------------------------|
| 1. Eyepiece | 8. Stage |
| 2. Rotary head | 9. Clips |
| 3. Objectives | 11. Condenser with iris diaphragm |
| 4. Estative | 12. Precondenser |
| 5. Focus limit screw | 13. Switch |
| 6. Coarse adjustment knob | 14. Light intensity adjustment knob |
| 7. Fine adjustment knob | |

3. TECHNICAL SPECIFICATIONS

Code	50116000
Head	Monocular
Eyepiece	WF10x / DIN 18mm
Nosepiece	Triple
Objectives	Achromatic: 4x, 10x, 40x (R)
Focusing	Coaxial coarse and fine focusing
Stage	Fixed with clamps (110x110 mm)
Condenser	Fixed with iris diaphragm and filter holder
Lighting	LED (5V, 1W)
Supply	Adapter 220V/50-60 Hz (included) or batteries 3xAA1.5V

These are the main features of the equipment:


- ◆ Metallic body very stable with knobs for coaxial coarse and fine focusing.
- ◆ Monocular head, inclined 45° and rotary 360°
 - Focusing cover: Independent 30mm with final stop
 - Fine 30 mm (0-200 micron each 2 microns)
- ◆ Fixed condenser with simple len (A.N._0.65) with iris diaphragm and mobile filter holder (blue filter Ø 32 mm included)
- ◆ Lighting high intensity LED type with switch and control potentiometer of the intensity.

Optical specifications

◆ Achromatic objectives: when talking about objectives, it is possible to distinguish between two types. Dry objectives are the ones where there is only air between the objective and the microscope slide, whereas immersion objectives are those that need immersion oil to be placed between the lens and the slide so as to obtain a higher luminosity.

Zuzi microscopes serie 116 are provided with three dry objectives (4x, 10x and 40x). The 40x objective is also retractile, which means that they present a built-in mechanism to avoid the breakage of the objective or the slide if the objective contacts the sample.

The features of objectives are inscribed on the barrel of each objective in the following way:

- 40:** Magnification
- 0.65:** Numerical aperture (N.A.)
- 160:** Mechanical tube length
- 0.17:** Cover glass thickness

Magnification	4x	10x	20x*	40x (R)
Working distance	37.50 mm	7.32 mm	1.26 mm	0.63 mm
Numerical aperture	0.10	0.25	0.40	0.65

(R) Retractable

* Optional

◆ Eyepieces: They are composed by two lenses separated by a diaphragm. Eyepieces work in combination with microscope objectives to further magnify the intermediate image so that sample details can be observed.

Type / magnification	WF/10x	WF/16x*
Focal distance	24.99 mm	15.58 mm
Field of view diameter	18	11

* Opcional

Attachable mechanical stage (code 90118590)

A mechanical stage can be attached on the stage of the Zuzi microscope model 116 to move the sample through the X and Y axes. The attachment present two lateral wheels and graduated axes to easily find structures of interest.

To put the attachable mechanical stage proceed as follows:

By using a screwdriver, remove the clips supplied with the fixed stage (**Figure 5**).



Figure 5

The attachable stage should fit in the holes of the stage marked as "b" (**Figure 6**). Put the attachable stage so as the screw of the attachment matches the small central hole at the upper part of the stage. The small metallic pieces at both sides of the screw must fit into the holes at both sides of the central hole of the stage (**Figure 6**).

Once the attachment is perfectly and securely placed on the stage, put the slide and hold it with the clamp of the attachment. By turning the lateral wheels of the attachment, the slide is moved over the stage.

Figura 6

similar to the objective's, thus avoiding marginal illumination that decreases contrast. Therefore, diaphragm must be more closed when using a lesser magnification objective with a small numerical aperture, and gradually open the diaphragm when user higher magnification objectives.

◆ The closing of the diaphragm by reducing the optical system's numerical aperture increases the focusing depth. Diffraction limit must not be surpassed as an excuse to increase the focusing depth.

◆ Diaphragm must not be used to reduce luminous diffraction. In spite of going straight on, light simple harmonic motion disperses when passing through a very closed diaphragm, provoking fuzzy vision and limiting the magnifying capacity of the microscope.

Accessories (not included)

1. Eyepiece video cameras

Code	Resolution	Output
59140052	640x480 pixels	Composite video
59140060	640x480 pixels	USB 2.0
59150060	1160x1024 pixels	USB 2.0
59150063	2048x1536 pixels	USB 2.0
59150065	2592x1944 pixels	USB 2.0
59150068	3264x2884 pixels	USB 2.0

2. Double head for demonstration (code 90119910)

This accessory allows two users to observe the sample at the same time.

3. Eyepieces and objectives

Eyepieces

Code	Description
90100126	Wide field eyepiece WF/16x

Achromatic objectives

Code	Description
90100203	20x

◆ Total magnification: It is the result of multiplying the eyepiece's magnification by the objective's magnification.

<i>Eyepieces</i>	<i>Objectives</i>	4x	10x	20x*	40x (R)
WF/10x	Total magnification	40	100	200	400
	Field of view	4.50 mm	1.80 mm	0.9 mm	0.45 mm
WF/16x*	Total magnification	64	160	320	640
	Field of view	2.75 mm	1.10 mm	0.55 mm	0.27 mm

(R) Retractable

* Optional

◆ Numerical aperture: It determines the objective's properties. The bigger it is, the brighter and better resolved the image will be.

◆ Working distance: it is the distance, in mm, between the microscope slide and the objective's front lens when the image is focused.

◆ Focal distance: It is the distance, in mm, from the main image plane to its image focus.

◆ Resolution: It is the reciprocal value of the separating power, that is to say, the smallest distance between two points on a sample that can still be distinguished as two separate entities.

◆ Field of view: Dimensions of the real field that is being observed, expressed in mm.

4. INSTALLATION / SETTING UP

Preliminary inspection

◆ Unwrap the microscope, take off the involving plastic and take off the polispán protection in which it comes fitted.

◆ Make sure that it does not present any damage because of the shipment. In case the equipment presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

◆ The returns have to be done before 15 days after the reception of the goods. Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- Microscope
- Achromatic objectives: 4x, 10x and 40x
- Eyepiece WF/10x
- Blue filter
- Plastic cover protector
- Adaptor
- Cleaner kit
- User's Manual
- Warranty Certificate

Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.


PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.

- ◆ Place the microscope on an horizontal, flat and stable table, making a free space for working.
- ◆ **VERY IMPORTANT:** Never hold the microscope by the stage nor the tube, as all the equipment's weight would fall on the fine adjustment screw and mechanical and precision parts would be slowly eroded.
- ◆ Avoid inflammable or toxic substances in the working area.
- ◆ The operator has to assume a comfortable, upright position with the back straight. It is convenient to work in a dark table in order to avoid an awkward fatigue and parasite light, as it may dazzle the operator or decrease the image sharpness.
- ◆ Please use the plastic cover provided whenever the microscope is not used so as to avoid dust laying on the optics.
- ◆ If you are not using the microscope for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

Setting up


◆ Lower the stage with the help of the coarse adjustment knob. Screw the objectives on the nosepiece so as when it is turned clockwise the objectives follow an ascendant way (4x-10x-40x) (Figure 2).

Figure 2

- ◆ The microscope is supplied with the eyepiece/s already placed. If changing the eyepieces is necessary, loosen the screw located at the ocular tube, change the eyepieces and fixed them again with the screw (Figure 3). The microscope head is normally placed facing the front of the microscope, but if it is necessary it can be placed in any other direction since it can be rotated in an angle of 360°.
- ◆ Connect the adaptor to supply bas. The unit works at 220 V, 50 Hz ± 10%.


Figure 3

Neither the manufacturer nor the distributor will assume any responsibility for damages caused in the equipment, installations or injuries if user does not follow electrical procedure. Electrical supply 220 V, 50 Hz ± 10%.

- ◆ Lower the stage according to the appropriate working distance for the objective with lesser magnification before placing the sample on it.
- ◆ Put the slide on the stage and, by looking directly at it (not by the eyepiece), place the sample centred on the stage opening.
- ◆ Looking through the eyepiece using the objective with lesser magnification, raise the stage by using the coarse adjustment knob until the image of the sample appears. Then, move the fine adjustment knob until the image is focused.
- ◆ Once the sample is focused with the objective with lesser magnification, screw the focus limit screw (Figure 4), so as avoiding the sample to hit accidentally the objectives and producing damages on them and on the sample.

Figure 4

◆ **IMPORTANT:** operation should always start by using the objective with lesser magnification, as it facilitates focusing and makes impossible to ruin the slides or get the objectives dirty. Furthermore, it is essential for adjusting the lighting previous to using the higher magnification, and it helps to see the whole topographic structure of the sample so as later concentrate on the more appealing parts to observe them at higher magnifications.

- ◆ By turning the nosepiece, pass to an objective with higher magnification. As they are parfocal objectives it will only be necessary to adjust focusing with the fine adjustment knob.
- ◆ Adjust both the diaphragm aperture until getting the light intensity desired. Normally, when operating with an intense lighting the condenser's numerical aperture should be lesser than the objective's, in order to avoid glares on the field of view that may distort the image contrast.
- ◆ Iris diaphragm enables decreasing the numerical aperture of the condenser to a value